

Implant/Abutment Interface Diameter

Our products are organized according to a “platform concept”, which facilitates treatment planning.

Narrow Platform (NP) for situations in which there is limited inter-dental space or a narrow crest.

Regular Platform (RP) for cases ranging from single anterior tooth loss to complete edentulism.

Wide Platform (WP) where additional loading can be expected or whenever a wider diameter implant/abutment is considered preferable.

6.0 Platform (6.0) where additional loading can be expected or whenever an even wider diameter implant/abutment is considered preferable.

External Connection

Brånemark System® Mk III Groovy*
 NobelSpeedy™ Groovy**
 Brånemark System® Mk III*
 Brånemark System® Mk IV
 Brånemark System® Zygoma TiUnite®
 Brånemark System® Zygoma

Internal Connection

NobelReplace™ Tapered Groovy
 NobelReplace™ Straight Groovy
 NobelSpeedy™ Replace
 Replace® Select Tapered
 Replace® Select Straight

One-Piece Implants

NobelDirect® Posterior***

*For RP Ø3.75 and RP Ø4 implants, use RP prosthetic components.

** For WP Ø5 and WP Ø6 implants, use WP prosthetic components.

***Use prosthetic components for Snappy Abutment™.

External Connection



NP Ø 3.5 mm interface



RP Ø 4.1 mm interface



WP Ø 5.1 mm interface



Internal Connection



NP Ø 3.5 mm interface



RP Ø 4.3 mm interface



WP Ø 5.0 mm interface



6.0 Ø 6.0 mm interface



NobelDirect® Posterior



RP Ø 4.3 mm



WP Ø 5.0 mm



6.0 Ø 6.0 mm

